PressCeramic

- [en] **Instructions for Use**
Lithium disilicate glass-ceramic ingot (LS2)
- [de] **Gebrauchsinformation**
Lithium-Disilikat-Glaskeramik-Rohling (LS2)
- [fr] **Mode d'emploi**
Lingotin vitrocéramique au disilicate de lithium (LS2)
- [it] **Istruzioni d'uso**
Grezzo in vetroceramica a base di disilicato di litio (LS2)
- [es] **Instrucciones de uso**
Lingote de cerámica vitrea de disilicato de litio (LS2)
- [pt] **Instruções de Uso**
Pastilha de vitrocerâmica de dissilicato de lítio (LS2)
- [sv] **Bruksanvisning**
Litiumdisilikat glaskeramisk puck (LS2)
- [da] **Bruksanvisning**
Lithium disilikat glaskeramik barre (LS2)
- [fi] **Käyttöohjeet**
Litiumdisilikaattilasikeraaminen nappi (LS2)
- [no] **Bruksanvisning**
Litiumdisilikat-glasskeram-rohling (LS2)
- [nl] **Gebruiksaanwijzing**
Lithiumdisilicaat keramische ingots (LS2)
- [el] **Οδηγίες Χρήσης**
Υαλοκεραμικός κύλινδρος διπυρρτικού λιθίου (LS2)
- [tr] **Kullanım Talimatları**
Lityum-disilikat cam seramik ingot (LS2)
- [ru] **Инструкция по применению**
Заготовка из дисиликата лития (LS2)
- [pl] **Instrukcja stosowania**
Pastylka z ceramiki szklanej dwukrzemowo litowej (LS2)
- [sl] **Navodila za uporabo**
Steklokeramični ingot iz litijevega disilikata (LS2)
- [hr] **Upute za uporabu**
Litij disilikatni staklokeramički valjčić (LS2)
- [cs] **Návod k použití**
Litij disilikátový sklokeramický ingot (LS2)
- [sk] **Návod na použitie**
Sklokeramický ingot z kremičitanu lítneho (LS2)
- [hu] **Használati útmutató**
Litium-diszilikát üveg-kerámia pellet (LS2)
- [sr] **Упутство за употребу**
Литијум-дисиликат стаклокерамички ингот (LS2)
- [mk] **Упатство за употреба**
Литиум дисиликат стаклокерамички ингот (LS2)
- [bg] **Инструкции за употреба**
Литиево дисиликатен стъклокерामичен слитък (LS2)
- [sq] **Udhëzimet e përdorimit**
Ingote qeramike qelqore Litium disilikat (LS2)
- [ro] **Instrucțiuni de utilizare**
Pastilă din ceramică vitroasă ranforsată cu disilicat de litiu (LS2)
- [uk] **Інструкція щодо використання**
Літію дісіліката склокерамічний інгот (LS2)
- [et] **Kasutamisujuhend**
Litium disilikaat klaaskeraamika valuplokk (LS2)
- [lv] **Lietošanas instrukcija**
Litija disilikāta stikla-keramikas lietnis (LS2)
- [lt] **Naudojimo instrukcija**
Ličio disilikato stiklo keramikos liejinys (LS2)
- [ar] **سببقة زجاج وخزف من ثنائي سيليكات الليثيوم (LS2)**
- [zh-tw] **二矽酸鋰玻璃陶瓷瓷塊 (LS2)**
- [zh-cn] **二硅酸鋰玻璃陶瓷瓷塊 (LS2)**
- [ja] **ニケイ酸リチウム・ガラスセラミックス・インゴット (LS2)**
- [ko] **리튬 디실리케이트 글라스-세라믹 잉곳 (LS2)**


Rx ONLY

CE 0123

Manufacturer:

Sagemax Bioceramics, Inc.
34210 9th Ave. South, Suite 218
Federal Way, WA 98003, USA
www.sagemax.com

DSSM AG
Im alten Riet 9
9494 Schaan/Liechtenstein

Date information
prepared:
2024-01-31 / Rev. 2

sagemax®

1 Intended use

Intended purpose

Single-tooth restorations in anterior and posterior teeth, 3-unit bridges up to the second premolar as the terminal abutment, implant-supported hybrid restorations for the replacement of single teeth.

Patient target group

Patients with permanent teeth, adult patients with dental implants

Intended users / Special training

- Dentists (clinical procedure), dental laboratory technicians (fabrication of restorations in the dental laboratory)
- No special training required.

Use

For dental use only.

Description

PressCeramic is a lithium disilicate glass-ceramic ingot (LS2) for the fabrication of fixed anterior and posterior restorations.

Technical data

Characteristics	Value
CTE (25 – 500 °C) [10 ⁻⁶ /K]	10.5 ± 0.5
Flexural strength (biaxial) [MPa]	≥ 360 / Typical average value: 470
Chemical solubility [µg/cm ²]	< 100
Type / Class	Type II / Class 3

According to ISO 6872:2015

Indications

Missing tooth structure in anterior and posterior teeth, partial edentulism in the anterior and posterior region

Types of restorations:

- Veneers
- Inlays
- Onlays (e.g. occlusal veneers, partial crowns)
- Crowns (on natural preparations and abutments)
- Abutment crowns
- Abutments
- Three-unit bridges up to the second premolar as the terminal abutment

Contraindications

- The use of the product is contraindicated if the patient is known to be allergic to any of its ingredients

Limitations of use

- Untreated bruxism (a splint is indicated after incorporation).
- Inlay bridges
- Cantilever bridges
- Adhesive bridges
- Restorations in the anterior region with a pontic width of > 11 mm
- Restorations in the premolar region with a pontic width of > 9 mm
- Temporary cementation of the restorations
- Very deep sub-gingival preparations
- Layering with a suitable veneering ceramic approved for the veneering of lithium disilicate glass-ceramic.
- A press furnace must be used that enables the operator to observe the indicated press parameters.

Additional limitations of use for minimally invasive anterior and posterior crowns

- Layer thicknesses below 1 mm
- Preparations with sharp edges
- Preparations that are not anatomically supported and feature varying layer thicknesses
- Conventional and self-adhesive cementation
- Build-up materials other than composite resin
- Absence of canine guidance
- Crowns on implants

Additional limitations of use for abutment restorations

-  Failure to observe the requirements of the implant manufacturer regarding the selected type of implant
- Failure to observe the manufacturer's instructions regarding the processing of the titanium bonding base
- Use of a luting composite other than Multilink Hybrid Abutment (Ivoclar Vivadent AG) for adhesive cementation
- Intraoral cementation of the ceramic structures to the titanium bonding base
- Temporary cementation of the crown on the hybrid abutment

Processing restrictions

-  Do not reuse.

Side effects

There are no known side effects to date.

Interactions

There are no known interactions to date.

Clinical benefit

Reconstruction of chewing function, restoration of esthetics

Composition

Dental ceramic material

After the manufacturing process of the glass-ceramic, a stable and inert network is formed in which the different elements are incorporated via oxygen bridges. The composition is determined as oxides.

Oxide	in wt.%
SiO ₂	57.0 – 80.0
Li ₂ O	11.0 – 19.0
K ₂ O	0.0 – 13.0
P ₂ O ₅	0.0 – 11.0
ZrO ₂	0.0 – 8.0
ZnO	0.0 – 8.0
other oxides and ceramic pigments	0.0 – 10.0

2 Application**Applications and processing techniques**

- Staining technique: Stain/Glaze firing using suitable staining and glazing materials.
- Cut-back technique: Incisal/vestibular veneering using a suitable veneering ceramic and subsequent Stain/Glaze firing.
- Layering technique: Dentin/Incisal firing using a suitable veneering ceramic and subsequent Stain/Glaze firing.

Ingot concept

PressCeramic ingots are offered as monochromatic ingots in four levels of translucency (MO, LT, MT, HT) and in special Opal shades in one size. For reasons of esthetics and the dental-lab protocol, the following processing techniques and types of restorations are recommended for the individual ingots depending on the degree of translucency.

	Translucency of the ingot				
	Opal	HT (High Translucency)	MT (Medium Translucency)	LT (Low Translucency)	MO (Medium Opacity)
Processing technique					
Staining technique	✓	✓	✓	✓	
Cut-back technique	✓	✓	✓	✓	
Layering technique					✓
Types of restorations					
Occlusal veneer ¹	✓	✓	✓		
Thin veneer ¹	✓	✓	✓		
Veneer	✓	✓	✓	✓	
Inlay		✓			
Onlay		✓	✓	✓	
Partial crown		✓	✓	✓	
Anterior and posterior crown			✓	✓	✓
3-unit bridge ²			✓	✓	✓
Hybrid abutment			✓	✓	✓
Hybrid abutment crown			✓	✓	

¹ The cut-back technique must not be used for thin and occlusal veneers.

² Only up to the second premolar as the distal abutment

Preparation instructions and minimum layer thicknesses

The preparation of the tooth structure is done according to the basic rules for all-ceramic restorations:

- No angles or edges
- Shoulder preparation with rounded inner edges and/or pronounced chamfer

When designing the restoration, the following minimum layer thicknesses (in mm) of the individual types of restorations and processing techniques must be observed:

Staining technique										
Types of restorations	Occlusal veneer	Thin veneer	Veneer	Inlay, onlay	Partial crown	Minimally invasive crown in the anterior and posterior region	Crown		Bridge	
							Anterior region	Posterior region	Anterior region	Premolar region
Incisal/occlusal	1.0	0.4	0.7	1.0 Fissure depth	1.0	1.0	1.5	1.5	1.5	1.5
Circular	1.0	0.3	0.6	1.0 Isthmus width	1.0	1.0	1.2	1.5	1.2	1.5
Connector dimensions	–	–	–	–	–	–	–	–	16 mm ² In general, the following applies: Height ≥ width	
Pontic width	–	–	–	–	–	–	–	–	≤ 11	≤ 9
Cementation	Mandatory adhesive cementation					Optional adhesive, self-adhesive or conventional cementation				

Cut-back technique										
Types of restorations	Occlusal veneer	Thin veneer	Veneer	Inlay, onlay	Partial crown	Crown		Bridge		
						Anterior region	Posterior region	Anterior region	Premolar region	
Incisal/occlusal	–	–	0.4	–	0.8	0.4	0.8	0.8	0.8	
Circular	–	–	0.6	–	1.5	1.2	1.5	1.2	1.5	
Connector dimensions	–	–	–	–	–	–	–	16 mm ² In general, the following applies: Height ≥ width		
Pontic width	–	–	–	–	–	–	–	≤ 11	≤ 9	
Cementation	Mandatory adhesive cementation					Optional adhesive, self-adhesive or conventional cementation				

Layering technique										
Types of restorations	Occlusal veneer	Thin veneer	Veneer	Inlay, onlay	Partial crown	Crown		Bridge		
						Anterior region	Posterior region	Anterior region	Premolar region	
Incisal/occlusal	–	–	–	–	–	0.6	0.8	0.8	0.8	
Circular	–	–	–	–	–	0.6	0.8	0.8	0.8	
Connector dimensions	–	–	–	–	–	–	–	16 mm ² In general, the following applies: Height ≥ width		
Pontic width	–	–	–	–	–	–	–	–	–	
Cementation	Mandatory adhesive cementation					Optional adhesive, self-adhesive or conventional cementation				

Dimensions in mm

Important


- The high-strength component (PressCeramic) of the restoration must always make up at least 50% of the total layer thickness of the restoration.
- In the fabrication of veneered or partially veneered restorations for large preparations, the excess space available must be filled by appropriately dimensioning the high-strength component (PressCeramic) and not by adding additional layering material.

The total layer thickness of the restoration is a combination of:

Total layer thickness of the restoration	0.8	1.0	1.2	1.5	1.8	2.0	2.5	3.0
Minimum thickness of the PressCeramic framework	0.4	0.5	0.6	0.8	1.0	1.1	1.3	1.6
Maximum layer thickness of the veneer	0.4	0.5	0.6	0.7	0.8	0.9	1.2	1.4


Dimensions in mm

Minimum layer thicknesses hybrid abutment

- The marginal shoulder width S_A must be at least 0.6 mm.
- Create an emergence profile with a right angle at the transition to the crown (see diagram).
- The wall thickness W_A must be at least 0.5 mm.
- The height H_A must not exceed twice the height of the titanium bonding base H_T .
- The hybrid abutment should be designed in a similar way as a prepared natural tooth:
 - Circular epi-/supragingival shoulder with rounded inner edges or a chamfer
 - In order for the crown to be cemented to the hybrid abutment using a conventional/self-adhesive cementation protocol, retentive surfaces and a sufficient "preparation height" must be observed.
- The width B_{AC} of the crown is limited to 6.0 mm from the axial height of contour to the screw channel of the hybrid abutment.

Minimum layer thicknesses hybrid abutment crown

- The marginal shoulder width S_A must be at least 0.6 mm.
- The wall thickness W_{AC} must be larger than 1.5 mm for the entire circumference.
- The opening of the screw channel must not be located in the contact point areas or areas with a masticatory function. If this is not possible, a hybrid abutment with a separate crown should be preferred.
- The width of the hybrid abutment crown B_{AC} is limited to 6.0 mm from the axial height of contour to the screw channel.
- The height H_{AC} must not exceed twice the height of the titanium bonding base by more than 2 mm.

Model and die preparation

Fabricate a model with detachable segments in the usual manner. Depending on the preparation, the spacer is applied to the die in several layers:

- For occlusal veneers, thin veneers, veneers, partial crowns, as well as single crowns, the spacer is applied in two layers up to max. 1 mm from the preparation margin (spacer application per layer 9–11 μ m).
- For inlays and onlays, the spacer is applied in up to three layers and up to the preparation margin.
- Also apply two layers for bridge reconstructions. Apply an additional layer at the intercoronal surfaces of the abutments (towards the pontic). This measure helps prevent undesired friction.
- For restorations on abutments, the procedure is the same as that on natural preparations.

Contouring

To contour the restoration, only organic waxes must be used, since they burn out without leaving residue. Contour the restoration in accordance with the desired processing technique (staining, cut-back, or layering technique). Please observe the following general notes for contouring:

- Observe the stipulated minimum layer thicknesses and connector dimensions of the respective indication and processing technique.
- Ensure precise contouring of the restorations, particularly in the area of the preparation margins. Do not over-contour the preparation margins, since this would require time-consuming and risky fitting procedures after pressing.
- For fully anatomical restorations, the possible occlusal relief must be taken into consideration as early as during the wax-up, since the application of the Stains and Glaze materials results in a slight increase in vertical dimensions.
- Do not model tips and edges with the cut-back and layering technique.
- For the layering technique, the frameworks should be anatomically reduced and designed to support the cusps.

Sprueing

Please observe the following notes when attaching the sprues to the wax-up:

- Before sprueing, weigh the ring base and record the weight.
- Always attach the sprues in the direction of flow of the ceramic and at the thickest part of the wax-up so that smooth flowing of the viscous ceramic during pressing is enabled.
- The attachment points must be rounded. Avoid angles and edges.
- Use a wax wire diameter between \varnothing 2.5 – 3 mm.
- Observe a wax wire length of min. 3 mm and max. 8 mm.
- Observe a distance of at least 3 mm between the objects.
- Observe a distance of at least 10 mm between the wax objects and silicone ring.
- The maximum length (wax objects + sprue) of 16 mm must not be exceeded.
- Attach the sprued wax-up at the "edge" of the investment ring base.
- Align the cervical margins of the wax objects to the silicone ring.

- To calculate the wax weight, weigh the loaded investment ring base again and then calculate the difference between the unloaded and loaded investment ring base.
- Use 1x 3 g ingot up to a maximum wax weight of 0.75 g.

Investing

-  Please follow the instructions for use of the investment material manufacturer.

Preheating

-  Please follow the instructions for use of the investment material manufacturer.
- Do not preheat the PressCeramic ingots and press plunger.

Pressing

-  Please follow the operating instructions of the press furnace manufacturer.

Basic procedure after completion of the preheating cycle:

- Remove the hot investment ring from the preheating furnace and place the cold PressCeramic ingot into it.
- Insert the ingots in the investment ring with the rounded, non-imprinted side facing down.
- Place the press plunger into the hot investment ring.
- Use the investment ring tongs to place the loaded investment ring upright and straight in the centre of the press furnace.
- Start press process with corresponding press parameters.

Press parameters

General press recommendations for PressCeramic

Ingot/ Translucency	Investment ring size	Start temperature	Heating rate	Pressing temperature	Holding time	Pressing time	Pressing pressure
	[g]	[°C/°F]	[°C/°F/min]	[°C/°F]	[min]	[min]	[N]
Opal, HT, MT	100	700/1292	60/108	920/1688	15	3	200 – 300 (approx. 3 – 4.5 bar/ 43 – 65 psi)
	200				25		
LT, MO	100	700/1292	60/108	925/1697	15	3	200 – 300 (approx. 3 – 4.5 bar/ 43 – 65 psi)
	200				25		

Procedure for determining the optimum pressing temperature

- Pin an inlay, a veneer and a crown and perform a test pressing.
- When all objects have been pressed out, gradually reduce the pressing temperature by 5 °C/9 °F until the pressing is incomplete.
- If not all objects have been pressed out, gradually increase the pressing temperature by 5 °C/9 °F until all objects have been pressed out.
- The lowest pressing temperature at which all objects are pressed out generally gives the best pressing results.

Important

- Depending on the press furnace used, the press temperature to be entered can sometimes deviate significantly from the recommended temperature. The recommended press temperatures should therefore only be seen as guidelines.
- Preheating furnace and pressing furnace must be calibrated regularly.
- The optimum pressing temperature depends on several factors. The use of reusable press plungers may require a 5 °C/9 °F higher pressing temperature. Depending on the investment material used, the press temperature can vary by ± 5 °C/± 9 °F. The higher the total liquid content of the investment material, the higher the press temperature is usually.

Divesting

Divest the investment ring as follows:

- Mark the length of the press plunger on the cooled investment ring.
- Separate the investment ring using a separating disc. This predetermined breaking point enables reliable separation of the press plunger and the ceramic material.
- Break the investment ring at the predetermined breaking point using a plaster knife.
- Always use polishing beads to divest the pressed objects (rough and fine divestment). Do not use Al₂O₃.
- Rough divestment is carried out with polishing beads at 4 bar (58 psi) pressure.
- Fine divestment is carried out with polishing beads at 2 bar (29 psi) pressure.
- Observe the blasting direction and distance to prevent damage to the object margins during divestment.

Finishing

Suitable grinding instruments are imperative for adjusting and finishing glass-ceramic materials. If unsuitable grinding instruments are used, chipping of the edges and local overheating may occur. The following procedure is recommended for finishing:

- Adjustment by grinding should be kept to a minimum.
- Overheating of the ceramic must be avoided. Low speed and light pressure must be observed.
- Separate the sprue using a suitable separating disc. Avoid overheating.
- Make sure that the minimum layer thickness of the restoration is maintained during finishing.
- Smooth out the attachment point of the sprue.
- Remove the spacer from the die. The restorations are tried in on the dies and carefully finished.
- Do not, under any circumstances, "post-separate" bridge framework connectors with separating discs. This may result in undesired predetermined breaking points, which will subsequently compromise the stability of the all-ceramic restoration.
- Check the occlusion and articulation and grind in the appropriate adjustments, if necessary.
- Design surface textures.
- Remove residues by briefly blasting the outer side of the restoration with Al₂O₃ 100 µm at 1 bar (14.5 psi) pressure and then steam-clean with a steam jet.

Completion

The press objects are finished using the staining, cut-back or layering technique with suitable ceramic materials. Please follow the instructions for use of the layering ceramic manufacturer.

Cementing

Preparation PressCeramic restoration	
Blasting	–
Etching	Etch the bonding surface for 20 sec with 5–9% hydrofluoric acid etching gel.  Please follow the instructions for use of the etching gel manufacturer.
Silanzing	Silanzize the bonding surface for 60 sec.  Please follow the instructions for use of the silane manufacturer.
Cementing	Mandatory adhesive cementation: occlusal veneer, thin veneer, veneer, inlay, onlay, partial crown Optional adhesive, self-adhesive or conventional cementation: crown, bridge  Please follow the instructions for use of the cementation manufacturer.

Bonding

	Preparation PressCeramic ceramic structure	Preparation titanium bonding base
Blasting	–	 Observe the manufacturer's directions
Etching	Etch the bonding surface for 20 sec with 5 – 9% hydrofluoric acid etching gel.  Please follow the instructions for use of the etching gel manufacturer.	–
Silanzing	Silanzize the bonding surface for 60 sec.  Please follow the instructions for use of the silane manufacturer.	Silanzize the bonding surface for 60 sec.  Please follow the instructions for use of the silane manufacturer.
Bonding	Multilink® Hybrid Abutment (Ivoclar Vivadent AG)  Please follow the instructions for use of the manufacturer.	

3 Cleaning, disinfection and sterilization

Hybrid abutments and hybrid abutment crowns must be cleaned, disinfected and sterilized immediately before use. Sagemax Bioceramics, Inc. recommends the following procedures:

Pre-cleaning

Sonicate the hybrid abutments and hybrid abutment crowns in water (minimum quality: drinking water) in an ultrasonic bath (e.g. Sonorex Digital 10P) for 2 min. Rinse under running tap water (minimum quality: drinking water) while brushing the inner and outer surfaces with a suitable brush (e.g. instrument cleaning brush with nylon bristles, Integra Miltex).

Cleaning and disinfection

Machine cleaning and disinfection in a washer-disinfector unit is preferred.

Machine cleaning and disinfection

If hybrid abutments and hybrid abutment crowns are exclusively machine cleaned, subsequent thermal disinfection is mandatory.

– Cleaning

Place the hybrid abutments and hybrid abutment crowns in a conventional sieve insert. Then place the sieve insert in a washer-disinfector compliant with ISO 15883 (e.g. Miele G7882, equipped with an upper basket Miele O 188/2). The automated cleaning procedure (e.g. based on Vario TD program) can be divided into the following steps:

- Cold water rinse for 5 min
- Cleaning at 50 ± 2 °C/122 ± 3.6 °F for 10 min with cleaning agent (e.g. neodisher MediZym, 0.2% v/v, Dr. Weigert)
- Cold water rinse for 2 min

– Disinfection

Thermal disinfection at 93 °C/199 °F for 5 min with demineralized water (A0 value > 3000 being reached at 90 °C/194 °F for 5 min). Only use deionized, low-germ (max. 10 cfu/ml) and low-endotoxin (max. 0.25 EU/ml) water for rinsing.

Manual cleaning and disinfection

– Cleaning

Immerse the hybrid abutments and hybrid abutment crowns in the cleaning agent (e.g. MD 520, undiluted) in an ultrasonic bath. Make sure that the immersed surfaces are completely covered in cleaning agent and that the cleaning agent is free of bubbles. After immersing the hybrid abutments and the hybrid abutment crowns in the ultrasonic bath, sonicate for 1 min.

Then clean the hybrid abutments and hybrid abutment crowns manually by thoroughly brushing the inner and outer surfaces with a suitable brush (e.g. instrument cleaning brush with nylon bristles, Integra Miltex) for at least 20 seconds per hybrid abutment and hybrid abutment crown until no residues are visible.

Subsequently rinse under running tap water (minimum quality: drinking water) for at least 10 seconds.

– Disinfection

Immerse the hybrid abutments and the hybrid abutment crowns in a disinfectant (e.g. MD 520, undiluted) in an ultrasonic bath and sonicate for 2 minutes. After sonication leave the hybrid abutments and the hybrid abutment crowns for 15 minutes at $20 \pm 2^\circ\text{C}/68 \pm 3.6^\circ\text{F}$ in the disinfectant. Make sure that the hybrid abutments and the hybrid abutment crowns are completely covered in disinfectant and that the disinfectant is free of bubbles.

Then immerse the hybrid abutments and the hybrid abutment crowns for 1 min in cold demineralized water to end the contact time with the disinfectant (this step is not a replacement for the extensive rinsing required for removing residues of the disinfectant when the hybrid abutments and the hybrid abutment crowns are conventionally cleaned).

After cleaning and disinfecting, thoroughly rinse the hybrid abutments and the hybrid abutment crowns with water. Only use deionized, low-germ (max. 10 cfu/ml) and low-endotoxin (max. 0.25 EU/ml) water for the final rinse.

Drying

Compressed air or clean, lint-free cellulose tissue.

Sterilization

Hybrid abutments and hybrid abutment crowns must be sterilized before use.

Sagemax Bioceramics, Inc. recommends one of the following sterilization procedures:

– For countries outside the United States:

To pack the items for sterilization, only use a sterile barrier system compliant with ISO 11607-1 (e.g. Steriking, Wipak) made of paper/film which is designated by the manufacturer for the extensive rinsing required for removing residues of the disinfectant when the hybrid abutments and the hybrid abutment crowns are conventionally cleaned. The sterile barrier system being used must be large enough. The filled sterile barrier system must not be stretched.

– For the United States:

Insert the products in a perforated basket with lid and wrap in two layers of 1-ply polypropylene wrap using sequential envelope folding technique prior to sterilization. Note: Users in the United States must make sure that the sterilizer and any sterilization accessories (e.g. sterilization wraps, pouches, basket, biological or chemical indicators) are approved by the FDA for the intended sterilization.

Steam-sterilize with a fractionated pre-vacuum process according to ISO 17665 in a steam sterilizer (e.g. Selectomat PL/666-1 CL) under the following conditions:

	Method	Conditions	Drying time
1	Steam sterilization (autoclave) Fractionated vacuum	134 °C (273 °F) for 4 min	Local practice
2	Steam sterilization (autoclave) [*] Fractionated vacuum	132 °C (270 °F) for 3 min	10 min
3	Steam sterilization (autoclave) ^{**} Fractionated vacuum	134 °C (273 °F) for 3 min	Local practice
4	Steam sterilization (autoclave) ^{***} Fractionated vacuum	134 °C (273 °F) for 18 min	Local practice

^{*} recommended for the USA

^{**} recommended for the UK

^{***} recommended for Switzerland and France

Storage

Sterilized products packed in a sterile barrier system (e.g. sterilization pouch) are intended for immediate use and must not be stored for longer than 48 hours.

4 Safety information

- In the case of serious incidents related to the product, please contact Sagemax Bioceramics, Inc., 34210 9th Ave. South, Suite 118, Federal Way, WA 98003, USA, website: www.sagemax.com, and your responsible competent authority.
- The current Instructions for Use are available in the download section of the Sagemax website (www.sagemax.com).
- Explanation of symbols: sagemax.com/eifu
- The Summary of Safety and Clinical Performance (SSCP) can be retrieved from the European Database on Medical Devices (EUDAMED) at <https://ec.europa.eu/tools/eudamed>. Basic UDI-DI: 084227112BABUT0077A

Warnings

- Observe the Safety Data Sheet (SDS, available on the website www.sagemax.com).
- The processing produces ceramic dust which may irritate the skin and eyes and which may result in lung damage. Make sure that the suction equipment at your workplace works flawlessly. Do not inhale grinding dust during finishing and wear a dust mask (particle class FFP2) as well as protective goggles.
- If applied frequently or for a prolonged period of time, peroxide-based professional bleaching agents (carbamide peroxide; hydrogen peroxide) as well as acidulated phosphate fluorides used for caries prevention can make the surface of existing PressCeramic restorations rough and matte.

Disposal information

Remaining stocks must be disposed of according to the corresponding national legal requirements.

Residual risks

Users should be aware that any dental intervention in the oral cavity involves a certain risks. Some of these risks are listed below:

- Chipping / fracture / decementation of the restorative material may lead to accidental swallowing or inhalation and dental re-treatment.
- Cement excess may lead to irritation of the soft tissue / gingiva. Progressive inflammation may lead to bone resorption or peri-implant disease.

5 Shelf life and storage

This product does not require any special storage conditions.

6 Additional information

Keep material out of the reach of children!

Not all products are available in all countries.

The product has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of use. The user is responsible for testing the product for its suitability and use for any purpose not explicitly stated in the Instructions.

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