NexxZr® Zirconia

NexxZr T, NexxZr T Multi, NexxZr+, NexxZr+ Multi, NexxZr+ Multi 2.0

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CE 0123

Manufacturer:

Sagemax Bioceramics, Inc.
34210 9th Ave. South, Suite 118
Federal Way, WA 98003, USA
www.sagemax.com

DSSM AG
Im alten Riet 9
9494 Schaan/Liechtenstein

Date information
prepared:
2026-04-24/Rev. 2

sagemax®

1 Intended use

Intended purpose

Crowns and bridges in the anterior and posterior region.

Patient target group

Patients with permanent teeth

Intended users / Special training

- Dentists (clinical procedure)
 - Dental laboratory technicians (fabrication of restorations in the dental laboratory)
- No special training required.

Use

For dental use only.

Description

NexxZr T, NexxZr T Multi, NexxZr+, NexxZr+ Multi, NexxZr+ Multi 2.0 is a yttrium-stabilized zirconium oxide for fixed all-ceramic dental restorations.

Note: The disc can be used to fabricate multiple restorations at different times.

Technical data

Material	NexxZr T	NexxZr T Multi	NexxZr+	NexxZr+ Multi	NexxZr+ Multi 2.0	
Shade range	White, OM1, OM2, OM3, A1, A2, A3, A3.5, A4, B1, B2, B3, B4, C1, C2, C3, C4, D2, D3, D4	BL1, BL2, BL3, BL4, A1, A2, A3, A3.5, A4, B1, B2, B3, B4, C1, C2, C3, C4, D2, D3, D4	White, A1, A2, A3, A3.5, B1, B2	A1, A2, A3, B1, B2, C2, D2	BL1, BL2, BL3, BL4, A1, A2, A3, A3.5, A4, B1, B2, B3, B4, C1, C2, C3, C4, D2, D3, D4	
Disc thicknesses	10, 12, 14, 16, 18, 20, 22, 25, 30 mm	16, 20, 25 mm	10, 12, 14, 16, 18, 20, 22, 25, 30 mm	16, 20 mm	14, 16, 20 mm	
Flexural strength [MPa]	Specifications	≥ 800	≥ 900	≥ 800	≥ 700	≥ 700
	Typical average value	1270	1200 (Dentin), 650 (Incisal)	1000 (white), 880 (pre-shaded)	850 (Dentin), 650 (Incisal)	850 (Dentin), 650 (Incisal)
Fracture toughness [MPa · m ^{1/2}]	>5.0	>5.0 (Dentin), >2.4 (Incisal)	>3.5	>3.5 (Dentin), >2.4 (Incisal)	>3.5 (Dentin), >2.4 (Incisal)	
CTE (25 – 500°C) [10 ⁻⁶ /K]	10.6 ± 0.5	10.5 ± 0.5	10.1 ± 0.5	10.4 ± 0.5	10.4 ± 0.5	
Chemical solubility [µg/cm ²]	< 100	< 100	< 100	< 100	< 100	
MDR classification	Class IIa	Class IIa	Class IIIa	Class IIa	Class IIa	
Type / Class ISO 6872:2015	Type II / Class 5	Type II / Class 5	Type II / Class 4	Type II / Class 4	Type II / Class 4	

Indications

- Missing tooth structure in anterior and posterior teeth
- Partial edentulism in the anterior and posterior region

Types of restorations:

Anterior and posterior restorations on prepared teeth and on approved implant abutment systems

Material	NexxZr T, NexxZr T Multi	NexxZr+, NexxZr+ Multi, NexxZr+ Multi 2.0
Full-contour crowns	✓	✓
Full-contour 3-unit bridges	✓	✓
Full-contour 4-unit and multi-unit bridges with max. 2 connected pontics	✓	–
Crown copings	✓	✓
3-unit bridge frameworks	✓	✓
4-unit and multi-unit bridge frameworks with max. 2 connected pontics	✓	–

Contraindications

- The use of the product is contraindicated if the patient is known to be allergic to any of its ingredients

Limitations of use

- Untreated bruxism (a splint is indicated after incorporation)
- The required minimum wall thicknesses and connector dimensions must be observed.
- Processing must be done with a milling system appropriate for the geometry of the disc.
- A sintering furnace must be used that enables the operator to observe the indicated sintering parameters.
- The final restoration must not be reused.

Side effects

There are no known side effects to date.

Interactions

There are no known interactions to date.

Clinical benefit

- Reconstruction of chewing function
- Restoration of esthetics

Composition

Material	NexxZr T	NexxZr T Multi	NexxZr+	NexxZr+ Multi	NexxZr+ Multi 2.0
Zirconium oxide (ZrO ₂)	> 87.0 %	> 87.0 %	> 85.0 %	> 85.0 %	> 85.0 %
Yttrium oxide (Y ₂ O ₃)	> 4.0 – ≤ 7.0 %	> 4.5 – ≤ 7.0 %	> 6.5 – ≤ 9.0 %	> 6.5 – ≤ 8.0 %	> 6.5 – ≤ 8.3 %
Hafnium oxide (HfO ₂)	≤ 5.0 %	≤ 5.0 %	≤ 5.0 %	≤ 5.0%	≤ 5.0%
Aluminium oxide (Al ₂ O ₃)	≤ 1.0 %	≤ 1.0 %	≤ 1.0 %	≤ 1.0%	≤ 1.0%
other oxides	≤ 1.0 %	≤ 1.5 %	≤ 1.0 %	≤ 1.0%	≤ 1.0%

2 Application

Application notes

I. Minimum thicknesses and connector dimensions

NexxZr T

Types of restorations	Anterior region		Posterior region	
	Minimum layer thickness in mm	Connector dimensions** in mm ²	Minimum layer thickness in mm	Connector dimensions** in mm ²
Crowns	0.4	–	0.6	–
3-unit bridges	0.6	7	0.6	9
4- and multi-unit bridges with 2 pontics***	0.6	9	0.7	12'
Cantilever bridges with one pontic	0.7	12'	0.7	12'

NexxZr T Multi

Types of restorations	Anterior region		Posterior region	
	Minimum layer thickness in mm	Connector dimensions** in mm ²	Minimum layer thickness in mm	Connector dimensions** in mm ²
Fully or partially anatomical restorations				
Crowns	0.8	–	1.0	–
3-unit bridges	1.0	9	1.0	9
4- and multi-unit bridges with 2 pontics***	1.0	9	1.0	12'
Cantilever bridges with one pontic	1.0	12'	1.0	12'
Frameworks; prerequisite: complete positioning in the dentin area				
Crowns	0.4	–	0.6	–
3-unit bridges	0.6	9	0.6	9
4- and multi-unit bridges with 2 pontics***	0.6	9	1.0	12'
Cantilever bridges with one pontic	1.0	12'	1.0	12'

NexxZr+, NexxZr+ Multi, NexxZr+ Multi 2.0

Types of restorations	Anterior region		Posterior region	
	Minimum layer thickness in mm	Connector dimensions** in mm ²	Minimum layer thickness in mm	Connector dimensions** in mm ²
Crowns	0.8	–	1.0	–
3-unit bridges	1.0	12'	1.0	16

* Height: 4 mm, width: 3 mm

** The specified minimum connector cross-section must be positioned in the dentin area of the disc

*** In Canada, bridge indications are limited to 6 units with a maximum of 2 connected pontics

II. Positioning guidelines

NexxZr T Multi

Disc height	16 mm	20 mm	25 mm
Incisal zone	3 mm	3 mm	3 mm
Transition zone	4 mm	4 mm	4 mm
Dentin zone	9 mm	13 mm	18 mm

NexxZr+ Multi

Disc height	16 mm	20 mm
Incisal zone	3.2 mm	4 mm
Transition zone	3.2 mm	4 mm
Dentin zone	9.6 mm	12 mm

NexxZr+ Multi 2.0

Disc height	16 mm	20 mm	25 mm
Incisal zone	3 mm	3 mm	3 mm
Transition zone	4 mm	4 mm	4 mm
Dentin zone	7 mm	9 mm	13 mm

General notes:

- The height positioning in the polychromatic disc determines the amount of incisal area in the restoration.
- Given the esthetic characteristics of single-tooth restorations, it is advisable to place the restoration about 1 mm below the top edge of the disc, in order to ensure a clearly visible incisal area.
- For fully or partially anatomical restorations, make sure that the specified minimum connector cross-section is positioned in the dentin area of the disc.
- Crown and bridge frameworks must be placed completely in the dentin area.

III. Nesting

Detailed information can be found in the Instructions for Use of the CAM software.

General notes:

- 3 holding bars must be attached to each single-tooth restoration.
- In multi-unit restorations, end units must have 2 holding bars attached (oral and vestibular). Attach holding bars to other units as required.
- The diameter of the holding bars must be at least 2.0 mm.
- The holding bars must be attached at least 1.0 mm above the preparation edge.
- The holding bars should be placed at the anatomical equator of the tooth so that no undercuts are created and the restoration can be easily processed from both sides.
- In the case of multi-unit restorations with a pronounced curvature, a sintering support structure is advisable.
- Provide end units with a holding bar vertical to the sintering support structure.
- Design the sintering support structure in an even thickness (2 – 5 mm).

IV. Milling

- For manual entry the specific shrinkage factor can be found on the disc. Detailed information can be found in the Instructions for Use of the CAM software.
- The imprinted side corresponds with the incisal / occlusal area.
- When securing the disc in the holder, care must be taken to ensure that the circumferential notch and the disc holder are absolutely clean and that the screws are tightened evenly crosswise.

V. Separating and finishing

- Grinding instruments which are approved for zirconium oxide (e.g. fine tungsten carbide burs) are recommended for separating the restorations. It is advisable to notch the holding bars on one side before the restoration is completely separated.
- Grinding instruments which are approved for zirconium oxide (e.g. fine tungsten carbide burs) are recommended to smooth out the attachment points of the holding bars. Rough tungsten carbide burs and/or grinding instruments are not suitable, as they may cause chipping, among other things.
- Remove any adhering zirconium oxide dust with a soft brush prior to sintering.

General notes:

- Non-sintered zirconium oxide restorations are susceptible to damage and fractures. This fact has to be kept in mind during the entire working procedure.
- Any adjustments should be carried out in the non-sintered state.
- In the non-sintered state, any contact with unsuitable liquids and liquids not approved for zirconium oxide (e.g. unpurified water and/or lubricant coolant) and/or contact media (e.g. occlusion spray) must be prevented.
- Use only light pressure for finishing.
- Do not, under any circumstances, post-separate bridge constructions with a separating disc. This may lead to predetermined breaking points in the area of the connectors and reduce the strength of the all-ceramic restoration.
- Rubber polishers should not be used as they condense the surface and cause contamination.
- Make sure that the minimum wall thicknesses of the restorations are maintained during finishing.
- Adhering zirconium oxide dust may be sintered to the restoration and lead to fitting inaccuracies.
- The non-sintered restoration must not be sonicated or steam-cleaned.
- The non-sintered restoration must not be blasted.

VI. Sintering

Sintering is a thermal process during which NexxZr T, NexxZr T Multi, NexxZr+, NexxZr+ Multi, NexxZr+ Multi 2.0 obtains its final physical and mechanical properties, such as its high strength and translucency.

In general, the following notes should be taken into consideration when sintering:

- Only fully dried restorations may be sintered in order to prevent damage to the furnace and/or the object.
- The use of sintering beads is not recommended.
- Sufficient atmospheric exchange in the sintering sagger must be ensured.
- The restorations must not come into contact with each other during sintering.
- The correct program selection must be observed.
- Sintering temperatures that are too low or too high and/or sintering times that are too short or too long will have a negative effect on the above mentioned final properties.
- Always keep the sintering accessories clean and free of dust so that the sintered restorations are not contaminated.

Information on sintering
NexxZr T, NexxZr+

	Number of units	Duration [h]	Phase	Temperature	Heating/ Cooling rate	Holding time
				[°C/°F]	[°C/°F/min]	[min]
Standard	1-5	~ 3.7	1	20-1300 / 68-2372	30 / 54	30
			2	1300-1530 / 2372-2786	40 / 72	60
			3	1530-900 / 2786-1652	15 / 27	-
			4	900-80 / 1652-176	20 / 36	-
	5-10	~ 5.2	1	20-1300 / 68-2372	30 / 54	60
			2	1300-1530 / 2372-2786	40 / 72	120
			3	1530-900 / 2786-1652	15 / 27	-
			4	900-80 / 1652-176	20 / 36	-
Long	1-20	~ 10.7	1	20-900 / 68-1652	10 / 18	10
			2	900-1530 / 1652-2786	3 / 5.4	150
			3	1530-80 / 2786-176	8 / 14.4	-
	> 21	~ 11.7	1	20-900 / 68-1652	10 / 8	10
			2	900-1530 / 1652-2786	3 / 5.4	210
			3	1530-80 / 2786-176	8 / 14.4	-
Over night	unlimited	~ 14.3	1	20-250 / 68-482	2 / 3.6	-
			2	250-1530 / 482-2786	4 / 7.2	240
			3	1530-80 / 2786-176	8 / 14.4	-

NexxZr T Multi, NexxZr+ Multi, NexxZr+ Multi 2.0

	Number of units	Duration [h]	Phase	Temperature	Heating/ Cooling rate	Holding time
				[°C/°F]	[°C/°F/min]	[min]
Speed	1-5	~ 4.8	1	20-1000 / 68-1832	60 / 108	10
			2	1000-1530 / 1832-2786	3 / 5.4	60
			3	1530-1100 / 2786-2012	50 / 90	-
			4	1100-100 / 2012-212	60 / 108	-
Long	unlimited	~ 9.6	1	20-900 / 68-1652	10 / 18	30
			2	900-1500 / 1652-2732	3.3 / 5.94	120
			3	1500-900 / 2732-1652	10 / 18	-
			4	900-300 / 1652-572	8.3 / 14.94	-

Processing after sintering

 Only use materials and tools which are approved for the processing and finishing of the zirconium oxide restorations. For detailed information, please refer to the Instructions for Use of the respective product.

The following aspects should be taken into consideration:

- Processing of sintered restorations should be kept to a minimum.
- The restoration should only be mechanically processed if it is absolutely necessary.
- Use only light pressure and low speed to adjust the restoration.
- Prevent sharp edges when finishing frameworks.
- Bridge connectors must not be post-separated with a separating disc.
- We recommend using a rubber polisher to smooth the basal side of the bridge connectors.
- The material-specific minimum wall thicknesses and connector dimensions must not be undermined during processing.
- Use only immaculate grinding instruments

- Remove zirconium oxide dust sintered to the restoration with suitable grinding instruments. As an alternative, zirconium oxide dust sintered to the restoration may be removed by blasting with Al₂O₃, 25 – 70 µm at a pressure of 1 bar (15 psi) or 70 – 110 µm at a pressure of 1.5 bar (22 psi).
- Rinse the restoration under running water or use the steam jet to remove any adhering residue and dry.
- Finish the restoration in accordance with the desired processing technique (staining, cut-back or layering).

Preparing for cementation

Condition the NexxZr T, NexxZr T Multi, NexxZr+, NexxZr+ Multi, NexxZr+ Multi 2.0 restoration by sandblasting the inner aspects of the crown with Al₂O₃, 25 – 70 µm, 1 bar (15 psi) or Al₂O₃, 70 – 110 µm, 1.5 bar (22 psi).

3 Safety information

- In the case of serious incidents related to the product, please contact Sagemax Bioceramics, Inc., 34210 9th Ave. South, Suite 118, Federal Way, WA 98003, USA, website: www.sagemax.com, and your responsible competent authority.
- The current Instructions for Use are available in the download section of the Sagemax website (www.sagemax.com).
- Explanation of symbols: sagemax.com/eifu
- The Summary of Safety and Clinical Performance (SSCP) can be retrieved from the European Database on Medical Devices (EUDAMED) at <https://ec.europa.eu/tools/eudamed>. Basic UDI-DI: 084227112ACERA0143P

Warnings

- Observe the Safety Data Sheet (SDS, available on the website www.sagemax.com).
- The processing of NexxZr T, NexxZr T Multi, NexxZr+, NexxZr+ Multi, NexxZr+ Multi 2.0 discs produces dust which may irritate the skin and eyes and which may result in lung damage. Make sure that the suction equipment at your milling machine and at your workplace works flawlessly. Do not inhale grinding dust during finishing and wear a dust mask (particle class FFP2) as well as protective goggles and gloves.

Disposal information

Remaining stocks must be disposed of according to the corresponding national legal requirements.

Residual risks

Users should be aware that any dental intervention in the oral cavity involves certain risks. Some of these risks are listed below:

- Chipping / fracture / decementation of the restorative material may lead to accidental swallowing or inhalation and dental re-treatment.
- Cement excess may lead to irritation of the soft tissue / gingiva. Progressive inflammation may lead to bone resorption or periodontitis.

4 Shelf life and storage

- In the original packaging
- In a dry place
- Do not expose to physical impact or vibrations.

5 Additional information

Keep material out of the reach of children!

Not all products are available in all countries.

The product has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of use. The user is responsible for testing the products for their suitability and use for any purpose not explicitly stated in the Instructions.

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